

### OPTIMUM RESOLUTION THANKS TO DIGITAL RAY-PATH® TECHNOLOGY

Digital Ray-Path® is an innovative calculation technique that uses a design engine to compensate the lens with a simulation of the binocular eye-lens system. Every unique lens is individually calculated guaranteeing an adapted solution for any prescription and base curve.



NON-PERSONALIZED PROGRESSIVE LENS



FULLY PERSONALIZED PROGRESSIVE LENS



### OPTIONS

#### MINIMUM FITTING HEIGHTS AVAILABLE

Ultimate N is available in 5 minimum fitting heights:

<b>MFH 14</b>	Minimum Fitting Height 14 mm
<b>MFH 15</b>	Minimum Fitting Height 15 mm
<b>MFH 16</b>	Minimum Fitting Height 16 mm
<b>MFH 17</b>	Minimum Fitting Height 17 mm
<b>MFH 18</b>	Minimum Fitting Height 18 mm

#### CUSTOMIZATION FOR CURVED FRAMES

For wrap frames, 100% personalization can be achieved by measuring the tilt of the lens.



A fully personalized progressive lens that offers an **extremely wide near visual field**




[www.iotamerica.com](http://www.iotamerica.com) / [www.digitalray-path.com](http://www.digitalray-path.com)  
3625 Del Amo Blvd., Suite 365, Torrance, CA 90503

Digital Ray-Path is a registered trademark of Indizen Optical Technologies.



Fully **personalized progressive** lens,  
enhanced for **reading**.



**PERSONALIZATION FOR EASY ADAPTATION AND VISUAL COMFORT**

Each lens is individually produced taking into consideration the parameters unique to every individual's face and frame combination. **Personalization is especially important for sport frames** to reduce the aberrations induced by the tilt position and curvature of the lens.



**PERSONALIZATION PARAMETERS**

It is essential to include **all personalization parameters** (described on the right page) **unique to each wearer's prescription data** when ordering an Ultimate N lens.



**A PREMIUM NEAR VISION PERSONALIZED PROGRESSIVE LENS**


The Ultimate N near visual field initiates a fast progression to provide wearers with a **generous zone for reading without having to constantly move their heads**. As a result, the near vision behavior of the Ultimate N lens is highly recommended for wearers that demand better near vision.




**ULTIMATE N DESIGN OVERVIEW**

The Ultimate N is a progressive lens developed to **improve vision when reading**, combined with wide, clear intermediate and distance visual fields. The main goal of this design is to enlarge the area devoted for near vision as much as possible, allowing the Ultimate N design to offer an **extremely wide reading zone**. The distance field of the lens behaves similar to a balanced design, making distance vision comfortable and easy to reach. Peripheral vision is sharp due to the design compensation based on the wearer's measured parameters.

**CYLINDER POWER MAP**



**DESIGN PERFORMANCE**



**PERSONALIZATION PARAMETERS**



**PRESCRIPTION & ADDITION**  
Digital Ray-Path<sup>®</sup> calculates and uses the accurate power that the user will truly perceive once the lenses are fitted on the frame.




**PANTOSCOPIC ANGLE**  
This is the angle in the vertical plane between the optical axis of a spectacle lens and the visual axis of the eye in primary position.



**NASOPUPILAR DISTANCE**  
The distance from the axis of symmetry of the face to the center of the pupil.



**WRAP ANGLE**  
The frame curvature.



**PUPILAR HEIGHTS**  
The vertical distance between the pupil center and the deepest part of the lens shape.



**BACK VERTEX DISTANCE**  
The distance between the cornea and the back surface of the lens.



**FRAME DIMENSIONS**  
Frame dimensions are used to calculate the final diameter, thickness of the lens and improve the efficiency of the optimization.



**NEAR WORKING DISTANCE**  
The distance from the lens to the typical reading position for the wearer.



Digital lens



Fully personalized



Digital Ray-Path<sup>®</sup>



Near vision enhanced



Multiple corridor



Short corridor available



Variable Inset



Wrap available